

**Racks for magnetic nanoparticle based separations****Bel-Art Products****Racks for magnetic nanoparticle based separations**

Rack suitable for all magnetic separation experiments involving affinity ligands separation of e.g. antibodies and nucleic acids.

- Model for 1.5 / 2.0 ml, 5 ml / 15 ml and 50 ml tubes
- Independent and rotatable sleeves with powerful magnets on 50 ml and 1.5/5 ml
- Easy removal of liquids
- 96-well PCR plate model has 24 magnets
- Flat bottom microplate model has magnets covering entire footprint to accept all standard micro-plates with any number of wells; 6, 12, 24, 48, 96 etc.
- Racks fit microplates following ANSI standards
- Clear acrylic allows unobstructed view of tube contents



Code	Description	Packaging
LLG06286784	Racks for magnetic nanoparticle based separations, For One 96-well PCR Tube Plate	1 pz.
LLG06286785	Racks for magnetic nanoparticle based separations, For One Standard Size Flat Microplate	1 pz.
LLG06286782	Racks for magnetic nanoparticle based separations, For Six 5 or 15ml Tubes	1 pz.
LLG06286783	Racks for magnetic nanoparticle based separations, For Ten 1.5 to 2.0ml Microcentrifuge Tubes	1 pz.
LLG06285493	Racks for magnetic nanoparticle based separations, For Three 50ml Tubes	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

