

**Magnetic stirrers BioStir**

DWK Life Sciences Inc.

Magnetic stirrers BioStir®

- Ideal for general mixing and difficult to dissolve solutes
- 150 to 1200 rpm, single and four place units
- Programmable, remote control and alarm capable

Specifications**Stirring quantity (H₂O): 3 l****Speed range: 150 to 1200 rpm****Permissible ambient temperature: 15 to 40 °C****Permissible relative moisture: 80 % to 31 °C / 50 % at 40 °C****Material: painted aluminium****Power supply: 200-240 VAC, 50/60 Hz****IP code: IP 21****Warranty: 2 years**

Code	Description	Packaging
LLG06286829	WHEATON-Bio-Stir®, 230VAC programmable magnetic stirrer single place, EU cable	1 pz.
LLG06286831	WHEATON-Bio-Stir®, 230VAC programmable magnetic stirrer with 4 places, EU cable	1 pz.
LLG06286832	WHEATON-Bio-Stir®, 230VAC programmable magnetic stirrer with 4 places, UK cable	1 pz.
LLG06286830	WHEATON-Bio-Stir®, 230VAC programmable magnetic stirrer single place, UK cable	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

