

**Dewar flasks, spherical, for LN2****KGW****Dewar flasks, spherical, for LN2**

DURAN®. DIN 16496. Spherical. Brushed aluminium casing with insulating lid and carrying handle.



Code	Description	Packaging
LLG09032115	Dewar flasks, spherical, for LN2, 5 Litres, Neck diam. 60 mm, Ext. diam. 260 mm	1 pz.
LLG09032741	Dewar flasks, spherical, for LN2, 3 Litres, Neck diam. 60 mm, Ext. diam. 225 mm	1 pz.
LLG09032740	Dewar flasks, spherical, for LN2, 1 Litres, Neck diam. 30 mm, Ext. diam. 175 mm	1 pz.
LLG09032127	Dewar flasks, spherical, for LN2, 10 Litres, Neck diam. 65 mm, Ext. diam. 330 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

