

**LLG-Dispenserbox, Acrylic Glass****LLG Labware****LLG-Dispenserbox, Acrylic Glass**

With hinged lid. For storage of laboratory tools, e.g. laboratory goggles, gloves or other small items.

Wall mounting possible (mounting material included)

-for about 12 goggles



Code	Description	Packaging
LLG06286593	Dispenserbox, Acrylic Glass, Dispenserbox, dim. (W x D x H) 216 x 216 x 200 mm	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

