

**Dewar vessels, cylindrical, for CO₂ and LN₂****KGW****Dewar vessels, cylindrical, for CO₂ and LN₂**

DURAN®. DIN 12492. Cylindrical. Borosilicate glass 3.3. With blue coated protective casing out of metal or structured aluminium casing. Large opening. With insulated lid. Carrying handles as indicated. For cooling and storage of temperature-sensitive goods e.g. with dry ice (CO₂). When storing liquid nitrogen (LN₂) a hole (diam. approx. 2mm) in centre of the lid is necessary.



Code	Description	Packaging
LLG09032131	Dewar vessels, cylindrical, for CO2 and LN2, 10 l, Int. diam. 200 mm, Int. height 350 mm, Figure a, With Handle	1 pz.
LLG09032132	Dewar vessels, cylindrical, for CO2 and LN2, 14 l, Int. diam. 200 mm, Int. height 500 mm, Figure a, With Handle	1 pz.
LLG09032742	Dewar vessels, cylindrical, for CO2 and LN2, 10 l, Int. diam. 200 mm, Int. height 350 mm, Figure d, With Handle	1 pz.
LLG09032133	Dewar vessels, cylindrical, for CO2 and LN2, 21 l, Int. diam. 250 mm, Int. height 480 mm, Figure c, With Side grips	1 pz.
LLG09032743	Dewar vessels, cylindrical, for CO2 and LN2, 21 l, Int. diam. 250 mm, Int. height 480 mm, Figure b, With Side grips	1 pz.
LLG06300165	Dewar vessels, cylindrical, for CO2 and LN2, 28 l, Int. diam. 250 mm, Int. height 620 mm, Figure b, With Side grips	
LLG09032134	Dewar vessels, cylindrical, for CO2 and LN2, 28* l, Int. diam. 250 mm, Int. height 620 mm, Figure c, With Side grips	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

