

**Extraction thimbles, Borosilicate glass****Hahnemühle****Extraction thimbles, Borosilicate glass**

Recommended for Soxhlet extraction with solvents, which are incompatible with cellulose fibres and for filtration of air and gas streams.

- Made of 100 % pure borosilicate fibres, without binding elements
- High loading capacity and highly permeable to air
- Thickness is 1.5 mm (diameter 33 mm)
- High retention of small particles, > 99 % according to BS 4400
- Good stability at high temperatures up to 500 °C
- Recommended for hot, moist, and acidic gases
- Gravimetric determination of particles in hot gases

Further sizes on request.



Code	Description	Packaging
LLG04006162	Extraction thimbles, Borosilicate glass, Ø 19 mm, Length 90 mm	25 pz.
LLG04006163	Extraction thimbles, Borosilicate glass, Ø 22 mm, Length 80 mm	25 pz.
LLG04006164	Extraction thimbles, Borosilicate glass, Ø 25 mm, Length 100 mm	25 pz.
LLG04006165	Extraction thimbles, Borosilicate glass, Ø 26 mm, Length 60 mm	25 pz.
LLG04006166	Extraction thimbles, Borosilicate glass, Ø 30 mm, Length 100 mm	25 pz.
LLG04006167	Extraction thimbles, Borosilicate glass, Ø 33 mm, Length 94 mm	25 pz.
LLG04006168	Extraction thimbles, Borosilicate glass, Ø 43 mm, Length 123 mm	25 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

