

**Removable Adapters (RN) for Luer Lock Hub Tubing****Hamilton****Removable Adapters (RN) for Luer Lock Hub Tubing**

For 250 µl RN needle syringes and larger. Port diameter 1.5 mm, Inert Gas Sampling Valve 1 mm.



Code	Description	Packaging
LLG06079139	Removable Adapters (RN) for Luer Lock Hub Tubing, Male Luer	1 pz
LLG06305865	Removable Adapters (RN) for Luer Lock Hub Tubing, Male Luer Lock	1 pz
LLG06306806	Removable Adapters (RN) for Luer Lock Hub Tubing, Female Luer	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

