

**Removable Adapters (RN) for Luer Lock Hub Tubing****Hamilton****Removable Adapters (RN) for Luer Lock Hub Tubing**

For 250 µl RN needle syringes and larger. Port diameter 1.5 mm, Inert Gas Sampling Valve 1 mm.



Code	Description	Packaging
LLG06306806	Removable Adapters (RN) for Luer Lock Hub Tubing, Female Luer	1 pz.
LLG06079139	Removable Adapters (RN) for Luer Lock Hub Tubing, Male Luer	1 pz.
LLG06305865	Removable Adapters (RN) for Luer Lock Hub Tubing, Male Luer Lock	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

