



Photometer photoLab and photoFlex, stored test protocols

Photometer photoLab® and photoFlex®, stored test protocols

A: Stored in photoLab® S6, S12, photoLab® 6000 Serie and photoLab® Spektral

B: Stored in photoLab® S12, photoLab® 6000 Serie and photoLab® Spektral

C: Stored in photoFlex® and photoFlex® Turb®

Reagent-free tests:

Cu - copper plating bath, 690nm, 820nm - copper plating bath - B

CrO3 - chromium plating bath - chromium plating bath - B

Ni - nickel plating bath - nickel plating bath - B

HZ - Hazen colour - Hazen colour - B

IFZ - iodine number, 340nm, 455nm - iodine number - B

m-1 (DFZ) - colouring - FB445 - B, C

m-1 (DFZ) - colouring - FB525 - A, C

m-1 (DFZ) - colouring - FB620 - B, C

FAU - turbidity 620 - T620 - B

E - extinction - extinction - A, B, C

SAK - Spectral absorption coefficient at 254 and 436nm photoLab® 6100 VIS (436nm) and UV-VIS photoLab® 6600 (254/436nm)

Further tests on demand, range information is based on max. area and can also be smaller for photoFlex®.



Code	Description	Packaging
LLG09920958	Test kit, Type A6/25, For A, C, NH ₄ -N, Ammonium, Unit 0.20 - 8.00 mg/l, No. of determinations 25	1 pz.
LLG09920945	Test kit Type 01744, For A, pH	1 pz.
LLG09920929	Test kit, Type 14752, For B, C, NH ₄ -N, Ammonium, Unit 0.01 - 3.00 mg/l, No. of determinations 500	1 pz.
LLG09920901	Test kit, Type 14566, For A, C, Zn, Zinc, Unit 0.20 - 5.00 mg/l, No. of determinations 25	1 pz.
LLG09920927	Test kit Type 14544, For A, C, NH ₄ -N, Ammonium	1 pz.
LLG09920928	Test kit, Type 14559, For A, NH ₄ -N, Ammonium, Unit 4.0 - 80.0 mg/l, No. of determinations 25	1 pz.
LLG09920888	Test kit, Type 14554, For A, C, Ni, Nickel, Unit 0.10 - 6.00 mg/l, No. of determinations 25	1 pz.
LLG09920930	Test kit, Type 14785, For B, C, Ni, Nickel, Unit 0.02 - 5.00 mg/l, No. of determinations 250	1 pz.
LLG09920931	Test kit, Type 14776/1, For B, C, NO ₂ -N, Nitrite, Unit 0.005 - 1.00 mg/l, No. of determinations 1000	1 pz.
LLG09920960	Test kit, Type N5/25, For A, C, NO ₂ -N, Nitrite, Unit 0.010 - 0.700 mg/l, No. of determinations 25	1 pz.
LLG09920925	Test kit, Type 14694, For A, O ₂ , Oxygen, Unit 0.5 - 12.0 mg/l, No. of determinations 25	1 pz.
LLG09920952	Test kit, Type 00607/1, For B, C, O ₃ , Ozone, Unit 0.010 - 4.0 mg/l, No. of determinations 200	1 pz.
LLG09920872	Test kit, Type 14833, For A, Pb, Lead, Unit 0.10 - 5.00 mg/l, No. of determinations 25	1 pz.
LLG09920941	Test kit, Type 14683, For A, RH, Residual hardness, Unit 0.50 - 5.00 mg/l Ca, No. of determinations 25	1 pz.
LLG09920900	Test kit Type 14394, For B, SO ₃ , Sulfite	1 pz.
LLG09920948	Test kit Type 01746, For B, SO ₃ , Sulfite	1 pz.
LLG09920898	Test kit, Type 14548, For A, C, SO ₄ , Sulfate, Unit 5 - 250 mg/l, No. of determinations 25	1 pz.
LLG09920899	Test kit, Type 14564, For A, SO ₄ , Sulfate, Unit 100 - 1000 mg/l, No. of determinations 25	1 pz.
LLG09920904	Test kit, Type 14697, For B, C, Ten-A, Surfactants anionic, Unit 0.05 - 2.00 mg/l, No. of determinations 25	1 pz.
LLG09920866	Test kit, Type 00861, For A, C, Zn, Zinc, Unit 0.025 - 1.000 mg/l, No. of determinations 25	1 pz.



Code	Description	Packaging
LLG09920902	Test kit, Type 14832, For B, Zn, Zinc, Unit 0.05 - 2.50 mg/l, No. of determinations 90	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

