



## Stainless steel filter manifold system

Sartorius Lab Instruments

### Stainless steel filter manifold system

Combisart® Multibranch System with modular structure: For microbiological analysis. Basic system without outlet for holding all dismountable Sartorius vacuum funnels either completely or by using appropriate adapters



Code	Description	Packaging
LLG09604216	Accessories for Filter Manifold system Combisart®, Base support, stainless steel for mounting a single funnel onto a manifold	1 pz.
LLG09604210	Accessories for Filter Manifold system Combisart®, Funnel, 100 ml	1 pz.
LLG09604211	Accessories for Filter Manifold system Combisart®, Funnel, 500 ml	1 pz.
LLG09604213	Accessories for Filter Manifold system Combisart®, Lid, stainless steel for 500 ml funnel	1 pz.
LLG09604215	Accessories for Filter Manifold system Combisart®, Silicone lid seal for 500 ml funnel	1 pz.
LLG09604227	Suction apparatusTriple	1 pz.
LLG09604228	Stainless steel filter manifold system	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

