



## Self-Seal Sterilization Pouches

### Self-Seal Sterilization Pouches

Self-seal sterilization pouches are pre-folded and assure accurate and fast closing without the requirement of a heat sealing device. They are especially suited for use in small general practices, dental surgeries and by users who want to avoid costs of additional sealing equipment. The self-seal sterilization pouches are constructed from transparent PET/PP multi-layer copolymer film, medical grade paper and adhesive strip which provides excellent resistance to over-exposure of steam and clean opening. Water based, non-toxic process indicators for steam and ethyleneoxide sterilization complying with ISO 11140-1 are applied on the paper surface and help to differentiate between processed and unprocessed packages.

- Reinforced film to avoid tearing during opening
- Clean and accurate indicator color change
- Tear, fibre-free opening and aseptic presentation
- Triple band seal for higher package integrity



Code	Description	Packaging
LLG06272283	Self-Seal Sterilization Pouches, dim. 130 x 290 mm	200 pz
LLG06272284	Self-Seal Sterilization Pouches, dim. 190 x 340 mm	200 pz
LLG06272285	Self-Seal Sterilization Pouches, dim. 300 x 430 mm	200 pz
LLG06272280	Self-Seal Sterilization Pouches, dim. 75 x 200 mm	200 pz
LLG06272281	Self-Seal Sterilization Pouches, dim. 90 x 230 mm	200 pz
LLG06272282	Self-Seal Sterilization Pouches, dim. 90 x 270 mm	200 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

