

**Shaker VORTEX 2****IKA****Shaker VORTEX 2**

Vortex shaker especially for short time operation (Touch function) or continuous operation.

Specifications

Type of movement: orbital

Orbital diameter: 4 mm

Max. shaken quantity: 400 g

Motor rating input / output: 39 / 9 W

Speed range: 500 to 2500 rpm

Dimensions (W x H x D): 120 x 138 x 140 mm

Weight: 3.9 kg

Enclosure according to DIN EN 60529: IP 21



Code	Description	Packaging
LLG07971443	Orbital shaker Vortex 2	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

