



Closures for narrow-mouth reagent bottles, series 310, UN-approved

Kautex

Closures for narrow-mouth reagent bottles, series 310, UN-approved

UN-approved.

Type 1: Light cap, white, PP, tamper-evident ring, white, PTFE cup-shaped liner (without figure)**Type 2: Light cap**, black, HDPE, tamper-evident ring, black, PTFE cup-shaped liner**Type 3: Bell cap**, black, PP, tamper-evident ring, red, cup-shaped PTFE liner**Type 4: Bell cap**, red, PP, with degassing membrane, tamper-evident ring, red, cup-shaped, perforated PTFE liner

Code	Description	Packaging
LLG09073380	Closures for narrow-mouth reagent bottles, series 310, UN-approved, Type 3, Bell cap, For bottles 500, 1000, 2500 ml, For Thread o.d. 45 mm	1 pz.
LLG09073381	Closures for narrow-mouth reagent bottles, series 310, UN-approved, Type 4, Bell cap, For bottles 500, 1000, 2500 ml, For Thread o.d. 45 mm	1 pz.
LLG09073379	Closures for narrow-mouth reagent bottles, series 310, UN-approved, Type 2, Light cap, For bottles 500, 1000, 2500 ml, For Thread o.d. 45 mm	1 pz.
LLG09073480	Closures for narrow-mouth reagent bottles, series 310, UN-approved, Type 1, Light cap, For bottles 100 ... 250 ml, For Thread o.d. 32 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

