



Undergloves white, polyester or nylon

Undergloves white, polyester or nylon

- The hem on the wrist is overlocked to prevent the edge from fraying
- Exceptionally breathable
- Stretchy and fit perfectly
- Colour Detection of the cuffs
- Carton with 30 x 10 pairs



Code	Description	Packaging
LLG04659018	Undergloves white, polyester or nylon, Size XL, Material polyester, Colour brown	300 pz.
LLG04659022	Undergloves white, polyester or nylon, Size XL, Material nylon, Colour brown	300 pz.
LLG04659020	Undergloves white, polyester or nylon, Size M, Material polyester, Colour green	300 pz.
LLG04659024	Undergloves white, polyester or nylon, Size M, Material nylon, Colour green	300 pz.
LLG04659019	Undergloves white, polyester or nylon, Size L, Material polyester, Colour grey	300 pz.
LLG04659023	Undergloves white, polyester or nylon, Size L, Material nylon, Colour grey	300 pz.
LLG04659021	Undergloves white, polyester or nylon, Size S, Material polyester, Colour white	300 pz.
LLG04659025	Undergloves white, polyester or nylon, Size S, Material nylon, Colour white	300 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

