



Cleanroom papers ASPURE

Cleanroom papers ASPURE

The cleanroom papers are available in 5 different colours.

- 72 g / m²
- Cellulose, water-soluble resin polymer
- Suitable for cleanroom class ISO 5
- Surface resistance: 1011W
- Packed dustfree in cleanroom



Code	Description	Packaging
LLG04659102	Cleanroom papers ASPURE, Size A4, Colour blue, Cleanroom Class ISO 5	2500 pz.
LLG04659103	Cleanroom papers ASPURE, Size A4, Colour white, Cleanroom Class ISO 5	2500 pz.
LLG04659104	Cleanroom papers ASPURE, Size A3, Colour blue, Cleanroom Class ISO 5	1250 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

