



Cleanroom Wipes ASPURE, polyester

Cleanroom Wipes ASPURE, polyester

Extremely low particle number and evaporation residues. Laser-cutting techniques minimise the formation of dust and lint on cutting edges. The use of polyester long fibers prevents the generation of lint.

- Knitted made of 100 % polyester
- Lasercut
- Cleaned with ultrapure water
- Double packed



Code	Description	Packaging
LLG04659070	Cleanroom Wipes ASPURE, polyester, dim. 229 x 229 mm, Cleanroom Class ISO 5, Package contents 1 x 150 pieces	150 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

