

**Sampler Mini-Sampler, PTFE****Bürkle****Sampler Mini-Sampler, PTFE**

For ultrapure sampling. Also suitable for viscous oils.

Complete with one PFA Sample Jar 180 ml, 5 m PTFE suction tubing (6 x 8 mm), one tubing cutter, one tubing weight (V2A/nickel-plated brass) and 20 seal-it safety seal labels in case.



Code	Description	Packaging
LLG09303975	Sampler Mini-Sampler, PTFE, Type Mini sampler	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

