



Cleanroom Wipes ASPURE, polyester / nylon

Cleanroom Wipes ASPURE, polyester / nylon

Split fibers ensure an good wiping ability.

- Knitted made of 80 % polyester / 20 % nylon

- Heat cut

- Double packed

- Weight: 210 ± 10 g/m²

- Thickness: 410 ± 30 μ m

- Water Absorbency: >412.32 g/m²

- Clean Pack: Cleanroom clean with ultrapure water



Code	Description	Packaging
LLG04659071	Cleanroom Wipes ASPURE, polyester / nylon, dim. 229 x 229 mm, Cleanroom Class ISO 5	100 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

