



Cleanroom Wipes ASPURE, polyester / nylon

Cleanroom Wipes ASPURE, polyester / nylon

Good dust suppression and wiping efficiency due to split fiber.

- Knitted made of 87 % polyester / 13 % nylon
- Lasercut
- Double packed
- Weight: 130 g/m²
- Thickness: 360 µm
- Water absorbency: > 195 g/m²
- Clean Pack: Cleanroom clean with ultrapure water



Code	Description	Packaging
LLG04659072	Cleanroom Wipes ASPURE, polyester / nylon, dim. 229 x 229 mm, Cleanroom Class ISO 5	150 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

