



## Cleanroom Wipes ASPURE, polyester / nylon

### Cleanroom Wipes ASPURE, polyester / nylon

Good dust suppression and wiping efficiency due to split fiber.

- Knitted made of 87 % polyester / 13 % nylon
- Lasercut
- Double packed
- Weight: 130 g/m<sup>2</sup>
- Thickness: 360 µm
- Water absorbency: > 195 g/m<sup>2</sup>
- Clean Pack: Cleanroom clean with ultrapure water



Code	Description	Packaging
LLG04659072	Cleanroom Wipes ASPURE, polyester / nylon, dim. 229 x 229 mm, Cleanroom Class ISO 5	150 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

