

**Cleanroom Wipes ASPURE, polyester / rayon****AS ONE****Cleanroom Wipes ASPURE, polyester / rayon**

Cleanroom wipes made of nonwoven fabric.

- 30 % polyester / 70 % rayon
- Quattro folded
- Weight: 30 g/m<sup>2</sup>
- Thickness: 360 µm
- Water absorbency: 226 g/m<sup>2</sup>



Code	Description	Packaging
LLG04659073	Cleanroom Wipes ASPURE, polyester / rayon, dim. 250 x 250 mm, Cleanroom Class ISO 6, Package contents 30 x 100 wipes	3000 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

