

**Beakers, PTFE****ISOLAB****Beakers, PTFE**

Manufactured from chemically inert PTFE with excellent chemical resistance, high mechanical strength and high temperature resistance between -200 to +260 °C. Flat base enables a flat contact with the surface for effective heat transfer.



Code	Description	Packaging
LLG06286607	Beakers, PTFE, Nominal capacity 500 ml, Ext. diam. 98 mm, Height 116 mm,, Spout +, Graduations -, Material PTFE, Type of beaker low form	1 pz.
LLG06286608	Beakers, PTFE, Nominal capacity 1000 ml, Ext. diam. 98 mm, Height 150 mm,, Spout +, Graduations -, Material PTFE, Type of beaker low form	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

