

**Beakers, stainless steel****ISOLAB****Beakers, stainless steel**

Manufactured from high quality, chemical and corrosion resistant AISI 304 grade stainless steel. Robust seamless construction guarantees a long life cycle. Highly polished smooth interior enables much easier cleaning.



Code	Description	Packaging
LLG06286609	Beakers, stainless steel, Nominal capacity 100 ml, Ext. diam. 50 mm, Height 60 mm,, Spout -, Graduations -, Material Stainless steel, Type of beaker low form	1 pz.
LLG06286612	Beakers, stainless steel, Nominal capacity 1000 ml	1 pz.
LLG06286610	Beakers, stainless steel, Nominal capacity 250 ml, Ext. diam. 65 mm, Height 85 mm,, Spout -, Graduations -, Material Stainless steel, Type of beaker low form	1 pz.
LLG06286611	Beakers, stainless steel, Nominal capacity 500 ml	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

