

**Measuring jugs with handle, stainless steel****ISOLAB****Measuring jugs with handle, stainless steel**

Manufactured from high quality, chemical and corrosion resistant AISI 304 grade stainless steel. Robust seamless construction guarantees a long life cycle. Highly polished smooth interior enables much easier cleaning. Moulded graduations on the side wall offer easy volume identification. Functional large spout minimize possible spills during liquid transfer.



Code	Description	Packaging
LLG06286613	Measuring jugs with handle, stainless steel, Nominal capacity 500 ml	1 pz.
LLG06286614	Measuring jugs with handle, stainless steel, Nominal capacity 1000 ml	1 pz.
LLG06286615	Measuring jugs with handle, stainless steel, Nominal capacity 2000 ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

