

**Screw caps, PP, GL 45****ISOLAB****Screw caps, PP, GL 45**

For perfect and tight closure of GL 45 screw necked laboratory bottles. Large ridges on the cap outer surface guarantees for increased grip and fast closure. Color coded caps enable fast, easy and practical identification of bottle content.



Code	Description	Packaging
LLG06286665	Screw caps, PP, GL 45, Colour Blue	10 pz.
LLG06286666	Screw caps, PP, GL 45, Colour Green	10 pz.
LLG06286663	Screw caps, PP, GL 45, Colour Red	10 pz.
LLG06286664	Screw caps, PP, GL 45, Colour Yellow	10 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

