

**LLG-Lab Chair****LLG Labware****LLG-Lab Chair**

Entry model for healthy sitting! The permanent contact backrest provides best support when working sitting down in the laboratory. The generously proportioned seat and backrest with integrated lumbar support provide additional support and therefore offer best laboratory ergonomics.

- Minimal joints: Hygienic design for easy cleaning and disinfection
- Resistant to disinfectants and chemicals: Resistant to all of the usual disinfectants and chemicals.
- Seat and backrest upholstery made of artificial leather (soft and comfortable, washable and antibacterial) or soft-touch PU foam (extremely robust, washable and antibacterial) in black.
- Quick adjustment: Every function can be operated comfortably and easy in sitting position.
- Seat height adjustment by gas spring.
- Permanent contact backrest.
- Backrest adjustable in height and inclination.
- Solid black plastic star base, flat design.
- The chrome-plated foot ring for the high chairs provides support through the wide tread area and can easily and quickly be adjusted in height.
- With certificate "Tested Safety" (GS mark).
- 10 years warranty, Made in Germany.

Accessories: Multifunctional armrests (pair): Adjustable in height, width and depth as well as swivelling.



Code	Description	Packaging
LLG06287753	Lab Chair, Type Lab Chair, Material Artificial Leather, Castors, Seat height adjustment 440-620 mm	1 pz.
LLG06287751	Lab Chair, Type Lab Chair, Material PU Foam, Glides, Foot Ring, Seat height adjustment 580-850 mm	1 pz.
LLG06287754	Lab Chair, Type Lab Chair, Material Artificial Leather, Glides, Foot Ring, Seat height adjustment 580-850 mm	1 pz.
LLG06287750	Lab Chair, Type Lab Chair, Material PU Foam, Castors, Seat height adjustment 440-620 mm	1 pz.
LLG06287752	Lab Chair, Type Lab Chair, Material PU Foam, Stop and Go Castors, Foot Ring, Seat height adjustment 620-890 mm	1 pz.
LLG06287755	Lab Chair, Type Lab Chair, Material Artificial Leather, Stop and Go Castors, Foot Ring, Seat height adjustment 620-890 mm	1 pz.
LLG06287756	Lab Chair, Type Armrests, Material,, Seat height adjustment - mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

