



Accessories for SCHOTT Cold light sources and EasyLED Spotlight

Accessories for SCHOTT Cold light sources and EasyLED Spotlight



Code	Description	Packaging
LLG09705221	Gooseneck light guide for SCHOTT Cold light sources, 1-branch, diam. 3.5 mm, L 500 mm, For KL 300 LED	1 pz.
LLG09705151	Gooseneck light guide for SCHOTT Cold light sources, 1-branch, diam. 4.5 mm, L 600 mm, For KL 1500 HAL, KL 1600 LED, KL 2500 LED	1 pz.
LLG09705168	Flexible light guide for SCHOTT Cold light sources, 1-branch, M-PVC coated, diam. 8.0 mm, L 1000 mm, For KL 1500 HAL, KL 1600 LED, KL 2500 LED	1 pz.
LLG09705220	Gooseneck light guide for SCHOTT Cold light sources, 2-branch, diam. 3.5 mm, L 500 mm, For KL 300 LED	1 pz.
LLG09705152	Gooseneck light guide for SCHOTT Cold light sources, 2-branch, diam. 4.5 mm, L 600 mm, For KL 1500 HAL, KL 1600 LED, KL 2500 LED	1 pz.
LLG09705226	Flexible light guide for SCHOTT Cold light sources, 2-branch, PVC coated, diam. 4.0 mm, L 1000 mm, For KL 300 LED	1 pz.
LLG06270979	Foot switch, For KL 2500 LED	1 pz.
LLG09705225	Flexible light guide for SCHOTT Cold light sources, 1-branch, PVC coated, diam. 4.5 mm, L 600 mm, For KL 300 LED	1 pz.
LLG07655134	Column adapter for base (Cat. No. 7.200 753), For KL 300 LED, EasyLED Spotlight	1 pz.
LLG09705177	Accessories for flexible and gooseneck light guides, Conversion filter, For Focussing attachment (Cat. No. 9.705 175)	1 pz.
LLG06301655	Day light filter, diam. 28 mm, For KL 1500 HAL	
LLG06270981	Diffusor, For EasyLED Ringlight (Cat. No. 6.258 176) / Ringlight Plus (Cat. No. 6.238 244) / VisiLED Ringlight (Cat. No. 6.231 031, 6.275 482)	1 pz.
LLG09705170	Accessories for flexible and gooseneck light guides, Focussing lens, with filter set, For Flexible / swan-neck light guide	1 pz.
LLG09705175	Accessories for flexible and gooseneck light guides, Focussing lens, without filter, For Flexible / swan-neck light guide	1 pz.
LLG09705196	Halogen bulb, 15 V / 150 W, For KL 1500 HAL	1 pz.
LLG09705272	Lens adapter, 58 mm to 66 mm, For All SCHOTT ringlights	1 pz.



Code	Description	Packaging
LLG06270980	Lens adapter, 60 mm to 66 mm, For All SCHOTT ringlights	1 pz.
LLG09705176	Accessories for flexible and gooseneck light guides, Polarisation filter, For Focussing attachment (Cat. No. 9.705 175)	1 pz.
LLG06270982	Polarisation filter set, For EasyLED Ringlight (Cat. No. 6.258 176) / Ringlight Plus (Cat. No. 6.238 244) / VisiLED Ringlight (Cat. No. 6.231 031, 6.275 482)	1 pz.
LLG06261529	Power cable for SCHOTT Cold light sources and Illumination systems, Power cord CH, 1.8 m	1 pz.
LLG06231176	Power cable for SCHOTT Cold light sources and Illumination systems, Power cord EU, 1.8 m	1 pz.
LLG06266145	Power cable for SCHOTT Cold light sources and Illumination systems, Power cord UK, 1.8 m	1 pz.
LLG06237972	Clamps for Cold light source KL 300 LED and EasyLED Spotlight, For Ø 32 mm column	1 pz.
LLG07910430	Clamps for Cold light source KL 300 LED and EasyLED Spotlight, For Ø 25 mm column	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
 Reagents operates with
 a Certified Quality
 Management System

