

**Jars with screw cap Nalgene, PC****Thermo Scientific****Jars with screw cap Nalgeneä, PC**

Cylindrical. Very durable. Stackable. Non-toxic. Suitable for cold storage rooms and refrigerators. Autoclavable - remove cap before autoclaving.

- Excellent for cold room and refrigerated use, easy to see contents
- Wide-mouth and straight sides provide full-access and complete sample retrieval ease


- 15, 30 and 60mL jars have natural polypropylene closures

- 125mL to 1L jars have white polypropylene closures with recessed tops for stacking

- 125 and 250mL sizes fit Osterizerä and other blenders that have Mason Jar threads, recommended for use with Fluid Transfer Closure 7.048 079

- Not liquid leakproof, not recommended for liquid transport

- Autoclavable/Transparent



Code	Description	Packaging
LLG09033239	Jars with screw cap Nalgene, PC, Nominal capacity 60 ml	4 pz.
LLG09033181	Jars with screw cap Nalgene, PC, Nominal capacity 15 ml	
LLG09033182	Jars with screw cap Nalgene, PC, Nominal capacity 30 ml	
LLG09033183	Jars with screw cap Nalgene, PC, Nominal capacity 60 ml	
LLG09033184	Jars with screw cap Nalgene, PC, Nominal capacity 125 ml	
LLG09033185	Jars with screw cap Nalgene, PC, Nominal capacity 250 ml	
LLG09033186	Jars with screw cap Nalgene, PC, Nominal capacity 500 ml	
LLG09033187	Jars with screw cap Nalgene, PC, Nominal capacity 1000 ml	1 pz
LLG09033237	Jars with screw cap Nalgene, PC, Nominal capacity 15 ml	4 pz.
LLG09033238	Jars with screw cap Nalgene, PC, Nominal capacity 30 ml	4 pz.
LLG09033240	Jars with screw cap Nalgene, PC, Nominal capacity 125 ml	4 pz.
LLG09033241	Jars with screw cap Nalgene, PC, Nominal capacity 250 ml	4 pz.
LLG09033242	Jars with screw cap Nalgene, PC, Nominal capacity 500 ml	4 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*



