



Refractometers, hand-held, precision

Refractometers, hand-held, precision

For simple, manual measurement of concentrations and refractive indices of liquids. Hand-held, precision refractometers with focussing eyepiece, soft rubber eyecups (to reduce straylight ingress, also spectacle lenses do not slip off and are not scratched), case and conversion tables. Hand-held refractometers cannot be fitted with a thermostatic controller and so measurements have to be temperature compensated.

Models with extra code M: Basic model

Models with extra code T: with automatic temperature compensation (ATC)

Models with extra code a: with automatic temperature compensation (ATC) and protection class IP65

Models with extra code H: for hot samples, with prism made of tempered glass



Code	Description	Packaging
LLG06230298	Refractometers, hand-held, precision, Type MASTER-20a, For Fruit juice, cutting oil	1 pz.
LLG06286461	Refractometers, hand-held, precision, Type MASTER-Honey/BX, For Honey	1 pz.
LLG09801301	Refractometers, hand-held, precision, Type MASTER-20T, For Brix	1 pz.
LLG06229977	Refractometers, hand-held, precision, Type MASTER-BAUME, For Wine	1 pz.
LLG09801372	Refractometers, hand-held, precision, Type MASTER-S10a, For Sodium chloride	1 pz.
LLG06231790	Refractometers, hand-held, precision, Type MASTER-TA, For Wine	1 pz.
LLG06286460	Refractometers, hand-held, precision, Type MASTER-KMW, For Wine	1 pz.
LLG06238823	Refractometers, hand-held, precision, Type MASTER-93H, For Brix	1 pz.
LLG07619590	Refractometers, hand-held, precision, Type MASTER 53a, For Brix	1 pz.
LLG09801472	Refractometers, hand-held, precision, Type MASTER-53T, For Brix	1 pz.
LLG09801359	Refractometers, hand-held, precision, Type MASTER-S/Mill-a, For Salinity, Sea water Spec. density	1 pz.
LLG06286459	Refractometers, hand-held, precision, Type MASTER-GOE, For Wine	1 pz.
LLG07970129	Refractometers, hand-held, precision, Type MASTER-OE, For Wine	1 pz.
LLG06235710	Refractometers, hand-held, precision, Type MASTER-20M, For Brix	1 pz.
LLG06261851	Refractometers, hand-held, precision, Type MASTER 500, For Universal application	1 pz.
LLG09801302	Refractometers, hand-held, precision, Type MASTER-53M, For Brix	1 pz.
LLG09801313	Refractometers, hand-held, precision, Type R-5000, For Refractive index	1 pz.
LLG09801333	Refractometers, hand-held, precision, Type HHR-2N, For Honey	1 pz.
LLG09801348	Refractometers, hand-held, precision, Type MASTER-BR, For Freezing point coolant*	1 pz.
LLG09801360	Refractometers, hand-held, precision, Type MASTER-RI, For Refractive index	1 pz.



Code	Description	Packaging
LLG09801370	Refractometers, hand-held, precision, Type MASTER-S10M, For Sodium chloride	1 pz.
LLG09801312	Refractometers, hand-held, precision, Type HSR-500, For Brix	1 pz.
LLG07652042	Refractometers, hand-held, precision, Type MASTER-100H, For Brix	1 pz.
LLG06230295	Refractometers, hand-held, precision, Type MASTER-80H, For Brix	1 pz.
LLG07653540	Refractometers, hand-held, precision, Type MASTER-S28a, For Sodium chloride	1 pz.
LLG09801335	Refractometers, hand-held, precision, Type MASTER-S28M, For Sodium chloride	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

