



Temperature probes, Pt100 for ama-digit ad 20 th

Temperature probes, Pt100 for ama-digit ad 20 th

Class A. Stainless steel sheathed probe, approx. 105 x 3mm L x dia.
With handle and approx. 1m cable.



Code	Description	Packaging
LLG09243114	Temperature probes, Pt100 for ama-digit ad 20 th, Type Pt100 air probe	1 pz.
LLG09243111	Temperature probes, Pt100 for ama-digit ad 20 th, Type Pt100 immersion probe	1 pz.
LLG09243112	Temperature probes, Pt100 for ama-digit ad 20 th, Type Pt100 penetration probe	1 pz.
LLG09243113	Temperature probes, Pt100 for ama-digit ad 20 th, Type Pt100 surface probe	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

