



## Flow cups with fixed nozzle

### Flow cups with fixed nozzle

With fixed nozzle made of stainless steel. Housing made of titanium anodised aluminum.



Code	Description	Packaging
LLG06226824	Flow cups with fixed nozzle, Flow nozzle 4 mm, DIN 53 211, Material Stainless Steel	1 pz.
LLG09149862	Flow cups with fixed nozzle, Flow nozzle 4 mm, DIN 53 211, Material Aluminium	1 pz.
LLG06088605	Flow cups with fixed nozzle, Flow nozzle 6 mm, ISO 2431, Material Aluminium	
LLG06901466	Flow cups with fixed nozzle, Flow nozzle 5 mm, ISO 2431, Material Aluminium	
LLG06091529	Flow cups with fixed nozzle, Flow nozzle 3 mm, similar to DIN 53 211, Material Aluminium	1 pz.
LLG06301841	Flow cups with fixed nozzle, Flow nozzle 6 mm, similar to DIN 53 211, Material Aluminium	1 pz.
LLG06302759	Flow cups with fixed nozzle, Flow nozzle 8 mm, similar to DIN 53 211, Material Aluminium	1 pz.



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

