



Flow cups with fixed nozzle

Flow cups with fixed nozzle

With fixed nozzle made of stainless steel. Housing made of titanium anodised aluminum.



Code	Description	Packaging
LLG06226824	Flow cups with fixed nozzle, Flow nozzle 4 mm, DIN 53 211, Material Stainless Steel	1 pz.
LLG09149862	Flow cups with fixed nozzle, Flow nozzle 4 mm, DIN 53 211, Material Aluminium	1 pz.
LLG06088605	Flow cups with fixed nozzle, Flow nozzle 6 mm, ISO 2431, Material Aluminium	
LLG06901466	Flow cups with fixed nozzle, Flow nozzle 5 mm, ISO 2431, Material Aluminium	
LLG06091529	Flow cups with fixed nozzle, Flow nozzle 3 mm, similar to DIN 53 211, Material Aluminium	1 pz.
LLG06301841	Flow cups with fixed nozzle, Flow nozzle 6 mm, similar to DIN 53 211, Material Aluminium	1 pz.
LLG06302759	Flow cups with fixed nozzle, Flow nozzle 8 mm, similar to DIN 53 211, Material Aluminium	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

