



## Tripod stand

### Tripod stand

Zinc-plated steel. With circular top and flared legs.

Zinc plating of tripod stands serves as temporary corrosion protection. May peel off after direct contact with open burner flame.



Code	Description	Packaging
LLG09034110	Tripod stand, Diam. 100 mm, Height 180 mm	1 pz
LLG06200910	Tripod stand, Diam. 120 mm, Height 250 mm	1 pz
LLG09034112	Tripod stand, Diam. 120 mm, Height 210 mm	1 pz
LLG09034114	Tripod stand, Diam. 140 mm, Height 220 mm	1 pz
LLG06059957	Tripod stand, Diam. 160 mm, Height 230 mm	1 pz



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

