



Safety Eyeshields uvex x-trend 9177

Safety Eyeshields uvex x-trend 9177

In addition to its stylish design, the uvex x-trend offers a high degree of safety with its large lens. It offers a clear view in all directions, 20% more cover than traditional spectacles, gives UV protection of up to 400 nm and comfortable temples.

uvex **supravision excellence** coated lenses are anti-fog (for at least 16 seconds) on the inside, while the outside offers extreme scratch and chemical resistance. The anti-fogging properties are permanent even after repeated cleaning. The lenses are also easy to clean and less susceptible to dirt due to the non-stick nanotechnology.



Code	Description	Packaging
LLG06280401	Safety Eyeshields uvex x-trend 9177, Colour black, Lens PC transparent / UV 2-1.2	
LLG09005149	Protection lenses x-trend 9177colour:	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

