



Stirrers for overhead stirrer EUROSTAR 20 high speed digital / control

IKA

Stirrers for overhead stirrer EUROSTAR 20 high speed digital / control



Code	Description	Packaging
LLG09197053	Stirrers for overhead stirrer EUROSTAR 20 high speed digital / control, Type R 1405, For Volume 0.25 - 30 l, Propeller, Rotor diam. 45 mm	1 pz.
LLG09197055	Stirrers for overhead stirrer EUROSTAR 20 high speed digital / control, Type R 1402, For Volume 1 - 30 l, Dissolver, Rotor diam. 42 mm	1 pz.
LLG09197054	Stirrers for overhead stirrer EUROSTAR 20 high speed digital / control, Type R 1401, For Volume 1 - 30 l, Propeller, Rotor diam. 55 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

