



Storage boxes, HDPE with tightly closing LDPE lid

Storage boxes, HDPE with tightly closing LDPE lid

- Natural
- Temperature range (without mechanical stress) from -20 to 95°C
- Acid and chemical resistant
- Suitable for cold storage rooms and refrigerators
- Stackable
- Food safe



Code	Description	Packaging
LLG09033270	Storage boxes, HDPE with tightly closing LDPE lid, 0.5 L, Width 103 mm, Depth 103 mm, Height 64 mm	1 pz.
LLG09033275	Storage boxes, HDPE with tightly closing LDPE lid, 1 L, Width 103 mm, Depth 103 mm, Height 127 mm	1 pz.
LLG09033280	Storage boxes, HDPE with tightly closing LDPE lid, 1 L, Width 103 mm, Depth 208 mm, Height 64 mm	1 pz.
LLG09033285	Storage boxes, HDPE with tightly closing LDPE lid, 1.5 L, Width 103 mm, Depth 208 mm, Height 94 mm	1 pz.
LLG09033286	Storage boxes, HDPE with tightly closing LDPE lid, 2 L, Width 103 mm, Depth 255 mm, Height 94 mm	1 pz.
LLG09033287	Storage boxes, HDPE with tightly closing LDPE lid, 2 L, Width 208 mm, Depth 208 mm, Height 64 mm	1 pz.
LLG09033288	Storage boxes, HDPE with tightly closing LDPE lid, 3.2 L, Width 208 mm, Depth 208 mm, Height 94 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

