

**Bottles, glass, culture medium, DURAN****DWK Life Sciences****Bottles, glass, culture medium, DURAN®**

Rimless. For glass caps. Glass Type I / neutral glass as per USP, EP and JP.
Autoclavable.



Code	Description	Packaging
LLG09010140	Bottles, glass, culture medium, DURAN®, 1000 ml, Diam. 105 mm, Neck diam. 46 mm, Height 237 mm	1 pz
LLG09010135	Bottles, glass, culture medium, DURAN®, 500 ml, Diam. 83 mm, Neck diam. 46 mm, Height 204 mm	1 pz
LLG09010130	Bottles, glass, culture medium, DURAN®, 300 ml, Diam. 71 mm, Neck diam. 31 mm, Height 169 mm	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

