



## Wire baskets, wire/nylon

### Wire baskets, wire/nylon

White. Nylon coated, strong wire.



| Code        | Description                             | Packaging |
|-------------|-----------------------------------------|-----------|
| LLG09033341 | Wire baskets,nylon-coated,100x100x100 m | 1 pz.     |
| LLG09033342 | Wire baskets,nylon-coated,120x120x120 m | 1 pz.     |
| LLG09033344 | Wire baskets,nylon-coated,140x140x140 m | 1 pz.     |
| LLG09033346 | Wire baskets,nylon-coated,160x160x160 m | 1 pz.     |
| LLG09033350 | Wire baskets,nylon-coated,200x200x200 m | 1 pz.     |



**CARLO ERBA Reagents S.A.S.**  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

