

**Clamps for Separating Funnels****LAUDA****Clamps for Separating Funnels**

Stainless steel, for screwing onto base trays 3950 and 3966.



Code	Description	Packaging
LLG09837970	Clamps for Separating Funnels, Type 3955, 250 ml, Max. clips per platform 4 / 8, Neck width 29 / 32 (ISO / Squibb / conical)	1 pz.
LLG09837975	Clamps for Separating Funnels, Type 3957, 50 ml, Max. clips per platform 6 / 11, Neck width 19 / 26 (ISO / Squibb)	1 pz.
LLG09837976	Clamps for Separating Funnels, Type 3958, 100 ml, Max. clips per platform 6 / 11, Neck width 19 / 26 (ISO / Squibb)	1 pz.
LLG09837977	Clamps for Separating Funnels, Type 3959, 250 ml, Max. clips per platform 4 / 8, Neck width 19 / 26 (ISO / Squibb)	1 pz.
LLG09837978	Clamps for Separating Funnels, Type 3956, 500 ml, Max. clips per platform 3 / 6, Neck width 29 / 32 (ISO / Squibb)	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

