

**Sterile stoppers SILICOSEN T-Type****Hirschmann****Sterile stoppers SILICOSEN® T-Type**

SILICOSEN® and BIO-SILICO® are specially processed culture plugs in sterile rubber with continuous bubbles for preparation, filling and sterilisation of culture media. Their uniform pore structure and good air permeability make them excellent for the cultivation of aerobic microorganisms.

The choice between SILICOSEN® and BIO-SILICO® culture stoppers depends on the requirements of the application. SILICOSEN® culture stoppers have low water evaporation which makes them particularly suitable for long term cultures. BIO-SILICO® culture stoppers have high air permeability, similar to cotton plugs. Both types are chemically resistant, temperature resistant, liquid repellent and reusable.

Moisture evaporation is minimal, thus there is less drying of the culture medium. Colour: natural.



Code	Description	Packaging
LLG09231017	Sterile stoppers SILICOSEN® T-Type, Type T-17, Cylinder length 10 mm, Cone length 30 mm, Joint size (NS) 13/16	1 pz.
LLG09231072	Sterile stoppers SILICOSEN® L-Type, Type L-12, Joint size (NS) 9/11, Cylinder length 40 mm, Cone length 10 mm	1 pz.
LLG09231010	Sterile stoppers SILICOSEN® T-Type, Type T-10, Cylinder length 13 mm, Cone length 17 mm, Joint size (NS) 6/9	1 pz.
LLG09231077	Sterile stoppers SILICOSEN® L-Type, Type L-17, Joint size (NS) 13/16, Cylinder length 30 mm, Cone length 20 mm	1 pz.
LLG09231015	Sterile stoppers SILICOSEN® T-Type, Type T-15, Cylinder length 15 mm, Cone length 20 mm, Joint size (NS) 11/14	1 pz.
LLG09231022	Sterile stoppers SILICOSEN® T-Type, Type T-22, Cylinder length 15 mm, Cone length 30 mm, Joint size (NS) 18/21	1 pz.
LLG09231019	Sterile stoppers SILICOSEN® T-Type, Type T-19, Cylinder length 12 mm, Cone length 30 mm, Joint size (NS) 15/18	1 pz.
LLG09231024	Sterile stoppers SILICOSEN® T-Type, Type T-24, Cylinder length 21 mm, Cone length 31 mm, Joint size (NS) 20/23	1 pz.
LLG09231028	Sterile stoppers SILICOSEN® T-Type, Type T-28, Cylinder length 20 mm, Cone length 35 mm, Joint size (NS) 20/23	1 pz.
LLG09231042	Sterile stoppers SILICOSEN® T-Type, Type T-42, Cylinder length 35 mm, Cone length 50 mm, Joint size (NS) 30/40	1 pz.
LLG09231052	Sterile stoppers SILICOSEN® T-Type, Type T-52, Cylinder length 40 mm, Cone length 60 mm, Joint size (NS) 40/50	1 pz.
LLG09231076	Sterile stoppers SILICOSEN® T-Type, Type T-76, Cylinder length 40 mm, Cone length 60 mm, Joint size (NS) 64/74	1 pz.
LLG09231038	Sterile stoppers SILICOSEN® T-Type, Type T-38, Cylinder length 30 mm, Cone length 35 mm, Joint size (NS) 26/35	1 pz.
LLG09231066	Sterile stoppers SILICOSEN® T-Type, Type T-66, Cylinder length 40 mm, Cone length 60 mm, Joint size (NS) 54/63	1 pz.



Code	Description	Packaging
LLG09231082	Sterile stoppers SILICOSEN® L-Type, Type L-22, Joint size (NS) 18/21, Cylinder length 30 mm, Cone length 30 mm	1 pz.
LLG09231012	Sterile stoppers SILICOSEN® T-Type, Type T-12, Cylinder length 14 mm, Cone length 17 mm, Joint size (NS) 9/11	1 pz.
LLG09231032	Sterile stoppers SILICOSEN® T-Type, Type T-32, Cylinder length 23 mm, Cone length 32 mm, Joint size (NS) 22/30	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

