



Swiss Boy laboratory jack accessories

Swiss Boy laboratory jack accessories

Swiss Boy laboratory stand accessories. Baseplates and ancilliary plates made of grey/black anodised aluminium. Support used to attach Swiss Boys in any position to apparatus scaffolding existing in many laboratories. Maximum support rod diameter 20mm. Black.



Code	Description	Packaging
LLG09117132	Swiss Boy laboratory jack accessories, Ancillary plate, Width 260 mm, Depth 180 mm	1 pz.
LLG09117130	Swiss Boy laboratory jack accessories, Ancillary plate, Width 200 mm, Depth 200 mm	1 pz.
LLG09117133	Swiss Boy laboratory jack accessories, Baseplate, Width 200 mm, Depth 200 mm	1 pz.
LLG09117131	Swiss Boy laboratory jack accessories, Ancillary plate, Width 300 mm, Depth 300 mm	1 pz.
LLG09117135	Swiss Boy laboratory jack accessories, Ancillary plate, Width 220 mm, Depth 320 mm	1 pz.
LLG09117134	Swiss Boy laboratory jack accessories, Ancillary plate, Width 500 mm, Depth 400 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

