


**Adhesive tape for label protection, PES and Tape dispenser, ABS****Bel-Art Products****Adhesive tape for label protection, PES and Tape dispenser, ABS****Adhesive tape for label protection**

Clear polyester tape with special synthetic adhesive. Resistant to moisture and most solvents. Test results prove the tape remains securely affixed with no discoloration after being autoclaved for 40 minutes at 125 °C, as well as after one week in liquid nitrogen at -196 °C. All tape sizes are 1 mil thick with 1 mil of acrylic adhesive. Excellent adhesion to glass, metal, PE etc. Supplied in rolls as outlined below.

Tape dispenser

Durable tape dispenser is molded of strong ABS plastic with a steel serrated cutting blade; tape roll held in place by aluminum roller and spring-loaded cover. Includes four non-slip rubber feet and screws. Dimensions (L x W x H): 206 x 117 x 60 mm.



Code	Description	Packaging
LLG09040550	Adhesive tape for label protection, PES and Tape dispenser, ABS, Tape dispenser, Width 117 mm, Length 206 m	1 pz.
LLG09040425	Adhesive tape for label protection, PES and Tape dispenser, ABS, Adhesive tape, Width 25 mm, Length 32,9 m	1 pz.
LLG09040438	Adhesive tape for label protection, PES and Tape dispenser, ABS, Adhesive tape, Width 38 mm, Length 32,9 m	1 pz.
LLG09040450	Adhesive tape for label protection, PES and Tape dispenser, ABS, Adhesive tape, Width 51 mm, Length 32,9 m	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

