

**Safety canisters for solvents****Rötzmeier****Safety canisters for solvents**

For safe storage, use and in-plant transports of flammable liquids. Explosion-proof by flame traps and pressure control valves in screw caps and metering devices (self-closing). From high quality and chemically resistant stainless steel 316/304, electro-polished, gaskets from PTFE, ergonomic handles. Space-saving shape, e.g. in safety cabinets.

Models with self-closing taps designed for horizontal storage and decanting.



Code	Description	Packaging
LLG09042911	Safety canisters for solvents, Type 10 KD, with self-closing, metering device 1 ½", 10 Litres, dim. (W x D x H) 130 x 280 x 380 mm, Weight 2.5 kg	1 pz.
LLG09042916	Safety canisters for solvents, Type 10 KK, with screw-cap 1 ½", 10 Litres, dim. (W x D x H) 130 x 280 x 380 mm, Weight 2.5 kg	1 pz.
LLG09042921	Safety canisters for solvents, Type 10 KZ*, with self-closing tap 1 ½", 10 Litres, dim. (W x D x H) 130 x 280 x 380 mm, Weight 2.5 kg	1 pz.
LLG09042910	Safety canisters for solvents, Type 05 KD, with self-closing, metering device 1 ½", 5 Litres, dim. (W x D x H) 130 x 200 x 310 mm, Weight 1.65 kg	1 pz.
LLG09042915	Safety canisters for solvents, Type 05 KK, with screw-cap 1 ½", 5 Litres, dim. (W x D x H) 130 x 200 x 310 mm, Weight 1.65 kg	1 pz.
LLG09042912	Safety canisters for solvents, Type 20 KDL*, with self-closing, metering device 1 ½", 20 Litres, dim. (W x D x H) 175 x 345 x 470 mm, Weight 3.7 kg	1 pz.
LLG09042917	Safety canisters for solvents, Type 20 KK, with screw-cap 1 ½", 20 Litres, dim. (W x D x H) 175 x 345 x 470 mm, Weight 3.7 kg	1 pz.
LLG09042941	Safety canisters for solvents, Type D 150, self-closing metering device 1 ½" with vent for 5 + 10 L, Litres, dim. (W x D x H) mm, Weight 0.17 kg	1 pz.
LLG09042943	Safety canisters for solvents, Type Z 150, self-closing tap 1 ½", Weight 0.17 kg	1 pz.
LLG09042920	Safety canisters for solvents, Type 05 KZ*, with self-closing tap 1 ½", 5 Litres, dim. (W x D x H) 130 x 200 x 310 mm, Weight 1.65 kg	1 pz.
LLG09042922	Safety canisters for solvents, Type 20 KZ*, with self-closing tap 1 ½", 20 Litres, dim. (W x D x H) 175 x 345 x 470 mm, Weight 3.7 kg	1 pz.
LLG09042976	Safety canisters for solvents, Type TR 2, Plug-in funnel, inclined, tube dia. 34mm, Litres, dim. (W x D x H) H: 120, Ø: 120 mm, Weight 0.35 kg	1 pz.



CARLO ERBA Reagents S.A.S.
Chaussée du Vexin – 27106 Val
de Reuil cedex
N° TVA: FR 63391048824
Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

