



Heating ovens, ED, FD, FED series - Avantgarde.Line

Heating ovens, ED, FD, FED series - Avantgarde.Line

The drying chambers with natural convection (Type designation ED) or forced convection (Type designation FD, FED) are renowned due to their new design, convenient operation, and efficiency. The latest APT.lineä technology offers outstanding temperature accuracy.

Equipment:

- Temperature range series ED: room temperature plus 5-6 °C to 300 °C
- Temperature range series FD, FED: room temperature plus 10 °C to 300 °C
- Controller with LCD display
- Electromechanical control of the exhaust air flap
- Class 2 independent temperature safety device (DIN 12880) with visual temperature alarm
- USB port for recording data

The number of shelves included in the scope of delivery is listed in the table.



Code	Description	Packaging
LLG06280069	Drying oven Model FD056, 60 ltr., USB port, 230 V 1 ~ 50/60 Hz	1 pz.
LLG06280295	Multifunctional oven FED 56, 56 L, 230V, 50/60Hz	1 pz.
LLG04658061	Drying ovens Model ED 56 up to +300°C, cap.57 ltrs 230 V, 50/60 Hz	1 pz.
LLG06280292	General purpose ovens/hot air sterilizers E 260, 260l, 230V, 50/60 Hz	1 pz.
LLG04658060	Drying ovens Model FD 115 up to +300°C, cap.116 ltrs 230 V 1~ 50/60 Hz	1 pz.
LLG06280284	Multifunctional heated/drying ovens FED 115, 115l,230 V.50/60 Hz.	1 pz.
LLG06280296	Multifunctional heated/drying ovens FED 260, 260l, 230 V., 50/60 Hz.	1 pz.
LLG06280362	Drying ovens Model FD 260, 260 ltr.,230 V 1N ~ 50/60 Hz	1 pz.
LLG04658059	Drying and heating chamber ED 115 with natural convection 114 l Interior volume, 230V	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

