

**Wide-mouth bottles, series 303, PP****Kautex****Wide-mouth bottles, series 303, PP**

Round, Natural without closure. With approximate graduations. Ideal for liquids, pastes and granulates. For taking samples, sending specimens and as transport containers. Can be filled with hot substances. Autoclavable to max. 121°C. **With screw closure.**

Note: Before autoclaving, unscrew closure and set lightly upside down on the container. Do not screw the closure down.



Code	Description	Packaging
LLG06287816	Wide-mouth bottles, series 303, PP, 1000 ml, Diam. 96 mm, Thread external diam. 65 mm, Height without closure 205 mm	1 pz.
LLG06287813	Wide-mouth bottles, series 303, PP, 100 ml, Diam. 48.5 mm, Thread external diam. 32 mm, Height without closure 94 mm	1 pz.
LLG06287814	Wide-mouth bottles, series 303, PP, 250 ml, Diam. 64.5 mm, Thread external diam. 40 mm, Height without closure 121 mm	1 pz.
LLG06287815	Wide-mouth bottles, series 303, PP, 500 ml, Diam. 77.5 mm, Thread external diam. 50 mm, Height without closure 155 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

