

**Dispensers, bottle-top, POLYFIX****Poulten & Graf****Dispensers, bottle-top, POLYFIX®**

For dispensing of all chemicals, where polypropylene is resistant. With glass dispensing barrel. Valve block not detachable for safety reasons.

-conformity-certified

-ergonomically designed plunger handle for comfortable operation

-user-friendly, volume setting with easy movement and rotatable fine adjustment

-PP valve block (GL32) with FIX adapter which aligns the dispenser on the bottle.

Supplied with 2 x PP bottle adapters (GL40 and GL45).



Code	Description	Packaging
LLG09283800	Dispensers, bottle-top, POLYFIX®, with glass piston and clear glass cylinder*, Volume 2 - 10 ml, Graduations 0.2 ml	1 pz
LLG09283801	Dispensers, bottle-top, POLYFIX®, with PTFE-coated piston and clear glass cylinder**, Volume 2 - 10 ml, Graduations 0.2 ml	1 pz
LLG09283802	Dispensers, bottle-top, POLYFIX®, with PTFE-coated piston and clear glass cylinder**, Volume 10 - 50 ml, Graduations 1.0 ml	1 pz
LLG09283803	Dispensers, bottle-top, POLYFIX®, with glass piston and amber dosing cylinder***, Volume 2 - 10 ml, Graduations 0.2 ml	1 pz
LLG09283804	Dispensers, bottle-top, POLYFIX®, with glass piston and amber dosing cylinder***, Volume 10 - 50 ml, Graduations 1.0 ml	1 pz
LLG09283810	Dispensers, bottle-top, POLYFIX®, with glass piston and clear glass cylinder*, Volume 10 - 50 ml, Graduations 1.0 ml	1 pz



CARLO ERBA Reagents S.A.S.  
 Chaussée du Vexin – 27106 Val  
 de Reuil cedex  
 N° TVA: FR 63391048824  
 Tél : +33 (0)2 32 09 20 00  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



CARLO ERBA  
 Reagents operates with  
 a Certified Quality  
 Management System