

**Aspirator carboys Nalgene, Type 2321, PP****Thermo Scientific****Aspirator carboys Nalgeneä, Type 2321, PP**

Aspirator carboy with screw cap and stopcock made of PP. With stainless steel handle. Graduated. Inside diameter of 86mm neck opening. Robust. Autoclavable only when empty, before autoclaving, unscrew cap.



Code	Description	Packaging
LLG09139921	Aspirator carboys Type 2321, PP, Type 2321, 9 Litres, Width 216 mm, Depth 146 mm, Height 356 mm	1 pz.
LLG09139922	Aspirator carboys Type 2321, PP, Type 2321, 20 Litres, Width 318 mm, Depth 219 mm, Height 381 mm	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

