

**Centrifuge bottles Nalgene, PC****Thermo Scientific****Centrifuge bottles Nalgeneä, PC**

PC with PP screw cap. Crystal clear. Good mechanical strength. Autoclavable. Leak-proof even at low speeds. In order to ensure that the system is leak-proof even at high speeds and when centrifuging dangerous substances, the 250ml and 1000ml sizes must be used in conjunction with sealing screw cap sets.



Code	Description	Packaging
LLG09315721	Centrifuge bottles Nalgene, PC, Volume 250 ml, Diam. 61.8 mm, Height 127.6 mm, Max. 27.500 rcf	1 pz.
LLG09315723	Centrifuge bottles Nalgene, PC, Volume 500 ml, Diam. 69.5 mm, Height 169.6 mm, Max. 13.700 rcf	1 pz.
LLG09315724	Centrifuge bottles Nalgene, PC, Volume 1000 ml, Diam. 97.6 mm, Height 188.4 mm, Max. 7.100 rcf	1 pz.
LLG09315725	Centrifuge bottles Nalgene, PC, Volume 1000 ml, Diam. 98.1 mm, Height 180.2 mm, Max. 7.100 rcf	1 pz



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

