

**Sleeve adapters, PTFE for Ultrasonic homogenisers SONOPULS****Bandelin****Sleeve adapters, PTFE for Ultrasonic homogenisers SONOPULS**

PTFE sleeve adapters for mounting probes into standard ground joint vessels.

For use with HD 2070.2 / 2200.2 / 4100 / 4200 / 4400 homogenisers.

NA 29 G: for vessels with standard ground joint 29/32 and for horns SH 70 / 213 G / 100 G / 200 G

NA 45 G: for vessels with standard ground joint 45/40 and for horns SH 70 / 213 / 219 / 225 G / 200 G / 400 G



Code	Description	Packaging
LLG09650228	Sleeve adapters, PTFE for Ultrasonic homogenisers SONOPULS, Type NA 29 G	1 pz.
LLG09650229	Sleeve adapters, PTFE for Ultrasonic homogenisers SONOPULS, Type NA 45 G	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

