



Glass sample vessels, Borosilicate glass 3.3 for Ultrasonic homogenisers SONOPULS

Glass sample vessels, Borosilicate glass 3.3 for Ultrasonic homogenisers SONOPULS

For direct ultrasonic irradiation.



Code	Description	Packaging
LLG09650235	Cooling vessel Type KG 3, For probes up to Ø 2 ... 13 mm, Volume 15 ml	1 pz.
LLG09650236	Type DG 3 flow-through vessel with cooling jacket , For probes up to Ø 2 ... 13 mm, Volume max. 5,6 l/h ml	1 pz.
LLG09650223	Cooling vessel KG 5, 70 ml	1 pz.
LLG09650224	Type DG 5 flow-through vessel with cooling jacket, For probes up to Ø 13 ... 25 mm, Volume max. 30 l/h ml	1 pz.
LLG09650268	Round anchor plate cell, glass RZ 2, 40ml	1 pz.



CARLO ERBA Reagents S.A.S.
 Chaussée du Vexin – 27106 Val
 de Reuil cedex
 N° TVA: FR 63391048824
 Tél : +33 (0)2 32 09 20 00
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

